

## Lamps

Bezeichnung	Modell	Art.-Nr.	OEM Art.-Nr.	Listpreis €
Deuterium Lamp	Holochrome	LD-GIL-100	HMD9510	€510,0
Deuterium Lamp	116/117/118/119	LD-GIL-101	100326	€350,0
Deuterium Lamp	151/152/155/156	LD-GIL-101	100326	€350,0
Deuterium Lamp	170 DAD	LD-GIL-105	2140-0590	€440,0
Xenon Lamp	122	LX-GIL-150MO	9021511	€425,0

## Pump Spares (All Models 301/303/305/306/307)

Bezeichnung	Type	Art.-Nr.	OEM Art.-Nr.	Listpreis €
Piston Seals (Black)	5SC	PS-GIL-102BK	5463088063	€36,3
	10SC & WSC	PS-GIL-202BK	5463125095	€36,3
	25SC & WSC	PS-GIL-302BK	5463050025	€46,2
	50SC	PS-GIL-402BK	5463070020	€50,6
	100SC	PS-GIL-502BK	5463394062	€40,7
	H1	PS-GIL-602BK	38013222	€60,5
	H2	PS-GIL-702BK	38013224	€60,5
	Piston Seals (Yellow)	5SC	PS-GIL-102YW	5463088863
10SC & WSC		PS-GIL-202YW	5463125895	€44,0
25SC & WSC		PS-GIL-302YW	5463050825	€60,5
50SC		PS-GIL-402YW	5463070820	€62,7
100SC		PS-GIL-502YW	5463394862	€37,4
H1		PS-GIL-602YW	38013220	€122,0
H2		PS-GIL-702YW	38013223	€60,5
Piston Assembly		5SC	PA-GIL-101	3650008
	10SC	PA-GIL-201	3650009	€113,3
	10WSC/WTI	PA-GIL-301	3650010	€136,4
	25SC	PA-GIL-401	3650011	€136,4
Inlet Check Valve Housing		CV-GIL-109	3645310	€73,7
Inlet Check Valve Cartridge	5/10/25SC & WSC	CV-GIL-103	36850082	€105,6
Inlet Check Valve Cartridge - Titanium (Long Life)	5/10/25SC & WSC	CV-GIL-103T		€101,2
Outlet Check Valve Housing		CV-GIL-110	3645309	€73,7
Outlet Check Valve Cartridge	5/10/25SC & WSC	CV-GIL-104	36850085	€105,6
Outlet Check Valve Cartridge - Titanium (Long Life)	5/10/25SC & WSC	CV-GIL-104T		€101,2
Outlet Washer	5/10/25SC & WSC	CW-GIL-113	3645329	€8,3
Outlet Frit Washer	5/10/25SC & WSC	CW-GIL-213	3645333	€12,7
Mixer Kit (new style)	811A	MK-GIL-116	AGT200S	€102,3
Replacement Mixer Kit Seal		MS-GIL-116		€8,3
Pump Return Spring		RS-GIL-116	36300058	€12,7
Anti Extrusion Gasket	5SC	AG-GIL-116	3650014K	€10,5
	10SC/WSC/WTI	AG-GIL-216	3650015K	€10,5
	25SC	AG-GIL-316	3650016K	€13,2

## For Solvent Lines

Bezeichnung	Art.-Nr.	OEM Art.-Nr.	Listpreis €
White Flangeless nuts 1/8" (pk10) Fits pump inlet & priming valve 1109	P-303X	P-303X	€13,0
Yellow Ferrules 1/8" for P303 nuts (pk10)	P-300X	P-300X	€13,0
PTFE tubing 1/8" OD (price per metre)	PT10	PT10	€3,3
Solvent Filter SS 10u 1/8" (inc. fitting)	A-302A	A-302A	€24,0
3 Way valve	001109	001109	€65,7
Female luer to fit 1109 valve	P-624	P-624	€4,0

**Tel: (02173) 8905 - 0** oder [sales@abimed.de](mailto:sales@abimed.de)

**online-shop [www.abimed.de](http://www.abimed.de)**

### DEC Caps

Bezeichnung	Art.-Nr.	OEM Art.-Nr.	Listpreis €
1ml Dec Cap for Aspec (PK/1000)	DC-GIL-001	2954698	€165,0
3ml Dec Cap for Aspec (PK/1000)	DC-GIL-003	2954699	€165,0

### After Pump Head (High Pressure End)

Bezeichnung	Art.-Nr.	OEM Art.-Nr.	Listpreis €
SS Fingertight 1/4"28 Thread (pk10) inc .ferrules Fits pumps/manometric module/mixer	F-160X	F-160X	€114,0
Replacement Peek ferrules for F160 nut (pk10)	F-162X	F-162X	€42,0
SS Tubing 0.020" ID (price per metre) Used from pumps to mixer	1T	1T	€8,0
SS Tubing 0.010" ID (price per metre) Used from mixer onwards	1U		€8,0
Peek tubing Orange 0.020" ID (price per metre) used from pumps to mixer	JR-610		€8,3
Peek tubing Blue 0.010" ID (price per metre) used from mixer onwards	JR-6009	JR-6009	€8,3
1 Piece extra long Peek Fingertights (pk10) Recommended	F-130X	F-130X	€54,0

#### Upchurch

Gerne stellen wir Ihnen unser komplettes Upchurch Sortiment vor. Rufen Sie uns an.

(02173) 8905 - 0

### Systems with Rheodyne Manual Injector 7125i and 7725i

Bezeichnung	Art.-Nr.	OEM Art.-Nr.	Listpreis €
Rhebuild Kit for 7125	7125-999		€151,0
Rhebuild Kit for 7725	7725-999		€151,0
5ul SS loop for 7125	7020		€41,0
10ul SS loop for 7125	7021		€40,0
20ul SS loop for 7125	7022		€36,0
50ul SS loop for 7125	7023		€38,0
100ul SS loop for 7125	7024		€40,0
2ul SS loop for 7725	7755-015		€112,0
5ul SS loop for 7725	7755-020		€41,0
10ul SS loop for 7725	7755-021		€38,0
20ul SS loop for 7725	7755-022		€36,0
50ul SS loop for 7725	7755-023		€38,0
100ul SS loop for 7725	7755-024		€40,0

#### Rheodyne

Gerne stellen wir Ihnen unser komplettes Rheodyne Sortiment vor. Rufen Sie uns an.

(02173) 8905 - 0

### Systems with Gilson Autosamplers

Bezeichnung	Modell	Art.-Nr.	OEM Art.-Nr.	Listpreis €
Short Septum Piercing Needle (requires sleeve)	231/232 A series	NP-GIL-201	27067131	€74,8
Short Needle Sleeve for NP-GIL-201		NSS-GIL-201	27067132	€94,6
Long Septum Piercing Needle (requires sleeve)	XL series	NP-GIL-202	27067235	€73,7
LongNeedle Sleeve for NP-GIL-202		NSL-GIL-202	27067236	€86,9
234 Needle - 1 piece		NP-GIL-101	27067251	€96,8
XL Needle - 1 piece		NP-GIL-102	27067238	€96,8
Injection Port Seal (old) (pk3)		IPS-GIL-001	2707243	€30,8
Injection Port Seal (new) (pk3)		IPS-GIL-002	27072755	€56,1
1. 1ml Coiled Transfer Tubing		CTT-GIL-001	2954583	€155,0
Aspec Non Septum Piercing Needle		NNA-GIL-101	27067374	€103,4
Aspec Septum Piercing Needle		NPA-GIL-102	27067383	€114,4
Aspec Injection Port Seal 1.5mm		IPS-GIL-003	2954674	€17,6
10ml Coiled Transfer Tubing for Aspec		CTT-GIL-010	2954672	€132,0
Calibrated Tubing SS with fittings		CTS-GIL-101	25051016	€12,3
Calibrated Tubing Peek		CTP-GIL-201	2707288	€68,2
500ul Syringe for 402 Dilutor		S402-GIL-005	25025347	€135,0
500ul Syringe for 401 Dilutor		S401-GIL-005		€208,0
1ml Syringe for 401 Dilutor		S401-GIL-010		€208,0